



Recombinant Escherichia coli Cysteine desulfurase (sufS)

Product Code	CSB-YP015769ENX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1LNI6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SMS-3-5 / SECEC)
Purity	>85% (SDS-PAGE)
Sequence	MKLPYLDYS ATTPVDPRVA EKMMQFMTMD GTFGNPASRS HRFGWQAEAA VDIARNQIAD LVGADPREIV FTSGATESDN LAIKGAANFY QKKGKHIITS KTEHKAVLDT CRQLEREGFE VTYLAPQRNG IIDLKELEAA MRDDTILVSI MHVNEIGVV QDIAAIGEMC RARGIYHVD ATQSVGKLPI DLSQLKVDLM SFGHKIYGP KGIGALYVRR KPRVRIEAQM HGGGHERGMR SGTLPVHQIV GMGEAYRIAK EEMATEMERL RGLRNRLWNG IKDIEEVYLN GDLEHGAPNI LNVSFNYVEG ESLIMALKDL AVSSGSACTS ASLEPSYVLR ALGLNDELAH SSIRFSLGRF TTEEEIDYTI ELVRKISIGRL RDLSPLEWEMY KQGVDLNSIE WAHH
Source	Yeast
Target Names	iscS
Protein Names	Recommended name: Cysteine desulfurase EC= 2.8.1.7
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.