



Recombinant Mycobacterium abscessus Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-EP018033MVF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1MH53
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium abscessus (strain ATCC 19977 / DSM 44196 / CIP 104536 / JCM 13569 / NCTC 13031 / TMC 1543)
Purity	>85% (SDS-PAGE)
Sequence	MARRPRPAGS EPDSGIPHIV ELVRSTIPPI HPAGLPFIAG GLGLAGLGVR NRWVRGTGLA LAGACAGFFR HPPRVPPNRA DVVVAADGQ VCLVDRAVPP PELGLPAEPL PRISIFLSVF DVHVQRVPVA GEATAVIHRS GQFLSADRAE ASVANERNSV QIRTRTGHDV IVVQIAGLIA RRIVCHAKVG DQLSIGDTYG LIRFG
Source	E.coli
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC=4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.