



# Recombinant Escherichia coli Queuine tRNA-ribosyltransferase (tgt)

<b>Product Code</b>	CSB-BP019146ENT
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B1XEZ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12 / DH10B)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKFELDTTDG RARRGRLVFD RGVVETPCFM PVGTYGTVKG MTPEEVEATG AQIILGNTFH LWLRPGQEIM KLHGDLHDFM QWKGPILTDS GGFQVFSLGD IRKITEQGVH FRNPINGDPI FLDPEKSMEI QYDLGSDIVM IFDECTPYPA DWDYAKRSME MSLRWAKRSR ERFDSLGNKN ALFGIIQGSV YEDLRDISVK GLVDIGFDGY AVGGLAVGEP KADMHRILEH VCPQIPADKP RYLMGVGKPE DLVEGVRRGI DMFDCVMPTR NARNGHLFVT DGVVKIRNAK YKSDTGPLDP ECDCYTCRNY SRAYLHHLDR CNEILGARLN TIHNLRYQR LMAGLRKAIE EGKLESFVTD FYQRQGREVP PLNVD
<b>Source</b>	Baculovirus
<b>Target Names</b>	tgt
<b>Protein Names</b>	Recommended name: Queuine tRNA-ribosyltransferase EC= 2.4.2.29 Alternative name(s): Guanine insertion enzyme tRNA-guanine transglycosylase
<b>Expression Region</b>	1-375
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.