



# Recombinant Escherichia coli ATP synthase subunit alpha (atpA)

<b>Product Code</b>	CSB-MP002344ENT
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B1X9W2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12 / DH10B)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQLNSTEISE LIKQRIAQFN VVSEAHNEGTV SVSDGVIR IHGLADCMQG EMISLPGNRY AIALNLERDS VGAVVMGPYA DLAEGMKVKC TGRILEVPVG RGLLGRVVNT LGAPIDGKGP LDHGDGFSAVE AIAPGVIERQ SVDQPVQTGY KAVDSMIPIG RGQRELIIGD RQTGKTALAI DAIINQRDSG IKCIYVAIGQ KASTISNVVR KLEEHGALAN TIVVVATASE SAALQYLAPY AGCAMGEYFR DRGEDALIIY DDLISKQAVAY RQISLLLRP PGREAFPGDV FYLHSRLLER AARVNAEYVE AFTKGEVKGK TGSLTALPII ETQAGDVSAF VPTNVISITD GQIFLETNLF NAGIRPAVNP GISVSRVGG AQT KIMK KLS GGIRTALAQY RELAAFSQFA SDLDDATRQK LDHGQKVTEL LKQKQYAPMS VAQQSLVLF A AERGYLADVE LSKIGSFEAA LLAYVDRDHA PLMQEINQTG GYNDEIEGKL KGILDSFKAT QSW
<b>Source</b>	Mammalian cell
<b>Target Names</b>	atpA
<b>Protein Names</b>	Recommended name: ATP synthase subunit alpha EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit alpha F-ATPase subunit alpha
<b>Expression Region</b>	1-513
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.