



Recombinant Human Selenoprotein W (SEPW1)

Product Code	CSB-YP021034HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P63302
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ALAVRVVYC GAUGYKSKYL QLKKKLEDEF PGRLDICGEG TPQATGFFEV MVAGKLIHSK KKG DG YVDTE SKFLKLVA AI KAALAAQ
Source	Yeast
Target Names	SELENOW
Protein Names	Recommended name: Selenoprotein W Short name= SelW
Expression Region	2-87
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a selenoprotein, which contains a selenocysteine (Sec) residue at its active site. The selenocysteine is encoded by the UGA codon that normally signals translation termination. The 3' UTR of selenoprotein genes have a common stem-loop structure, the sec insertion sequence (SECIS), that is necessary for the recognition of UGA as a Sec codon rather than as a stop signal. This protein shows highest expression in skeletal muscle and heart, and may be involved in oxidation-reduction reactions. A retroprocessed pseudogene, SEPW1P, has been identified and mapped to chromosome 1p35-34.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.