



Recombinant Callithrix jacchus Lys-63-specific deubiquitinase BRCC36 (BRCC3)

Product Code	CSB-EP002798CYL-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B0KWU8
Product Type	Recombinant Protein
Immunogen Species	Callithrix jacchus (White-tufted-ear marmoset)
Purity	>85% (SDS-PAGE)
Sequence	AVQVVQAVQ AVHLESDAFL VCLNHALSTE KEEVMGLCIG ELNDDTRSDS KFAYTGTEMR TVAEKVDAVR IVHIHSVIL RRS DKRKDRV EISPEQLSAA STEAERLAE L TGRPMRVVGW YHSHPHITVW PSHVDVRTQA MYQMMDQGFV GLIFSCFIED KNTKTGRVLY TCFQSIQAQK SSELHGPRD FWSSSKHISI EGQKEEERYE RIEIPIHVP HVTIGKVCLE SAVELPKILC QEEQDAYRRI HSLTHLDSVT KIHNGSVFTK NLCSQMSAVS GPLLQWLEDR LEQNQQHLRE LQQEKEELMQ ELSSLE
Source	E.coli
Target Names	BRCC3
Protein Names	Recommended name: Lys-63-specific deubiquitinase BRCC36 EC= 3.4.19.- Alternative name(s): BRCA1-A complex subunit BRCC36 BRCA1/BRCA2- containing complex subunit 3 BRCA1/BRCA2-containing complex subunit 36 BRISC complex subunit BR
Expression Region	2-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.