



Recombinant Escherichia coli Nicotinate phosphoribosyltransferase (pncB)

Product Code	CSB-MP015451ENX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B1LJT2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SMS-3-5 / SECEC)
Purity	>85% (SDS-PAGE)
Sequence	MTQFASPVLH SLLDTDAYKL HMQQAVFHYY YDVHVAAEFR CRGDDLLGIY ADAIREQVQA MQHLRLQDDE YQWLSALPFF QADYLNWLRE FRFNPEQVTV SNDNGKLDIR LSGPWREVIL WEVPLLAVIS EMVHRYRSPQ TDVAQALDTL ESKLVDFSAL TAGLDMSRFH LMDFGTRRRF SREVQETIVK RLQQESWFGV TSNYDLARRL SLTPMGTQAH EWFQAHQQIS PDLANSQRAA LAAWLEEYPD QLGIALTDCI TMDAFLRDFG VEFASRYQGL RHDSGDPVEW GEKAI AHYEK LGIDPQSKTL VFSDNLDLRK AVELYRHFSS RVQLSFGIGT RLTC DIPQVK PLNIVIKLVE CNGKPVAKLS DSPGKTICH D KAFVRALRKA FDLPHIKKAS
Source	Mammalian cell
Target Names	pncB
Protein Names	Recommended name: Nicotinate phosphoribosyltransferase Short name= NAPRTase EC= 2.4.2.11
Expression Region	1-400
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.