



Recombinant *Paeonia mascula* subsp. *russoi* Maturase K (matK)

Product Code	CSB-BP300137PAK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68754
Product Type	Recombinant Protein
Immunogen Species	<i>Paeonia mascula</i> subsp. <i>russoi</i> (Peony) (<i>Paeonia russoi</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEKSQGYLEL DKSWRHNFLY PLIFQEYIYA LAHEQGLNRS ILLENTDHDN KYSSLIVKRL ITLIHQQNHF LIFDNDNSQN PFWKHNNNLY SQTISEGFVI IVEIPFSPRF VDSLEEKKKI LKSNLRSIH SIFPFLEDQI LHLNFVANIL IPYPIHLEIV VQSLRYRVKD ASSLHLLRFF LFTLNKSISS FSKRNQRFFL FLYNHSHVY EY ESTFLFLRNK TSHLRSTSSG AFLERIFFYG KIKHLIEVFA NDFQAILWLF KDPFMHYVRY QGKSILASKR TSLRMNKWKY YLVNFWQCQF YVWSQPGRVS INQLSNHSLD FLGYLSSVRR NPLAVRSQML ENSFLTDNAI KKFDIIVLLI SLIGSLAKAK FCNVLGHPLS KPARADSSDS DIERFVRIC RNLSHYHSGS SKKKSlyRIK YILRLSCART LARKHKTTVR SFLKRLGSEL LEEFLTEDGQ VISLIFPRTS STSWRLYRGG IWYLDITCIN DLANHE
Source	Baculovirus
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.