



Recombinant Hepatitis delta virus genotype II Small delta antigen

Product Code	CSB-EP300229HFW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P69619
Product Type	Recombinant Protein
Immunogen Species	Hepatitis delta virus genotype II (isolate Japanese S-1) (HDV)
Purity	>85% (SDS-PAGE)
Sequence	MSQSETRRGR RGTREETLEK WITARKKAAEE LEKDLRKTRK TIKKLEENP WLGNIIVGIIR KGKDGEGAPP AKRPRTDQME VDSGPGKRPH KSGFTDKERE DHRRRKALEN KKKQLSAGGK ILSKEEEEEEL RRLTDEDEER KRRVAGPRVG DVNPSRGGPR GAPGGGFVPQ MAGVPESPFS RTGEGLDIRG TQGFP
Source	E.coli
Protein Names	Recommended name: Small delta antigen Short name= S-HDAg Alternative name(s): p24
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.