



# Recombinant Mouse Complexin-1 (Cplx1)

<b>Product Code</b>	CSB-MP005893MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P63040
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEFVMKQALG GATKDMGKML GGDEEKDPDA AKKEEERQEA LRQAEERKA KYAKMEAERE VMRQGIRDKY GIKKKEEREA EAQAAMEANS EGSLTRPKKA IPPGCGDEPE EEDESILDTV IKYLPGPLQD MFKK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Cplx1
<b>Protein Names</b>	Recommended name: Complexin-1 Alternative name(s): 921-S Complexin I Short name= CPX I Synaphin-2
<b>Expression Region</b>	1-134
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Proteins encoded by the complexin/synaphin gene family are cytosolic proteins that function in synaptic vesicle exocytosis. These proteins bind syntaxin, part of the SNAP receptor. The protein product of this gene binds to the SNAP receptor complex and disrupts it, allowing transmitter release.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.