



# Recombinant Mycoplasma pneumoniae Putative type-1 restriction enzyme specificity protein MPN\_089 (MPN\_089)

<b>Product Code</b>	CSB-EP300503MLW
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P75604
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGRIKTYDFD GEYVTWTTTRW SYAGSIYYRN GKFSASSNCG ILKVLNKEIN PKFLAYALKK EAKKFNNTTS AIPILRTQKV VEIPIDFPPL QIQEKIATIL DTFTELSAEL SAELSAELSA ELSAELRERK KQYAFYRDYL LNLKNWKEEN KYYKLGEIAQ KVLVGGEKPA DFSKEKNEVY KYPILSNNSK AEEFLVYSKT FRVEEKSITV SARGTIGAVF YRDFAYLPAV SLICFVPKEE FDIRFLFHAL RAIKFKKQGS ATGQLTVAQF KEYGIHVPSL KKQKEIAAIL DPLYSFFTDL NEGIPAEIEL RKKQLDYYQN FLFNWVQNQK AASIL
<b>Source</b>	E.coli
<b>Target Names</b>	MPN_089
<b>Protein Names</b>	Recommended name: Putative type-1 restriction enzyme specificity protein MPN_089 Alternative name(s): S.MpnORFAP Type I restriction enzyme specificity protein MPN_089 Short name= S protein
<b>Expression Region</b>	1-335
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.