



Recombinant Escherichia coli Uncharacterized protein ybiU (ybiU)

Product Code	CSB-BP300511ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75791
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MASTFTSDTL PADHKAAIRQ MKHALRAQLG DVQQIFNQLS DDIATRVAEI NALKAQGDV WPVLSYADIK AGHVTAEQRE QIKRRGCAVI KGHFPREQAL GWDQSM LDYL DRNRFDEVYK GPGDNFFGTL SASRPEIYPI YWSQAQMQR QSEEMANAQS FLNRLWTFES DGKQWFNP DV SVIYPDRIRR RPPGTTSKGL GAHTDSGALE RWLLPAYQ RV FANVFNGNLA QYDPWHAHR TEVEEYTVDN TTKCSVFRTF QGWTALSDML PGQGLLHVVP IPEAMAYVLL RPLDDVPED ELCGVAPGRV LPVSEQWHPL LIEALTSIPK LEAGDSVWWH CDVIHSVAPV ENQQGWGNVM YIPAAPMCEK NLAYAHKVKA ALEKGASPGD FPREDYETNW EGRFTLADLN IHGKRALGMD V
Source	Baculovirus
Target Names	ybiU
Protein Names	Recommended name: Uncharacterized protein ybiU
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.