



Recombinant Escherichia coli Uncharacterized protein ymfN (ymfN)

Product Code	CSB-EP300525ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75978
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSTKLTGYVW DGCAASGMKL SSVAIMARLA DFSNDEGVCW PSIETIARQI GAGMSTVRTA IARLEAEGWL TRKARRQGDG SSPHCAVVDE YEHATDALY TTMLTGMGAR RQPLMWAIT AGYNIEGPCY DKRREVIEML NGSVPNDEL GIIYTVDEGD DWTDPQVLEK ANPNIGVSVY REFLLSQQQR AKNNARLANV FKTKHLNIWA SARSAYFNLV SWQSCEDKSL TLEQFEGQPC ILAFDLARKL DMNSMARLYT REIDGKTHYY SVAPRFWVPY DTVYSVEKNE DRRTAERFQK WVEMGVLTVT DGAEVDYRYI LEEAKAANKI SPVSESPIDP FGATGLSHDL ADEDLNPVTI IQNYTNMSDP MKELEAAIES GRFHHDGNPI MTWCIGNVVG KTIPGNDDVV KPVKEQAENK IDGAVALIMA VGRAMLYEKE DTLSDHIESY GIRSL
Source	E.coli
Target Names	ymfN
Protein Names	Recommended name: Uncharacterized protein ymfN
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.