



Recombinant Escherichia coli Uncharacterized protein ynbC (ynbC)

Product Code	CSB-EP300535ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76092
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MENSRIPGEH FFTTSDNTAL FYRHWPALQP GAKKVIVLFH RGHEHSGRLQ HLVDELAMPD TAFYAWDARG HGKSSGPRGY SPLSARSVRD VDEFVRF AAS DSQVGLLEEVV VIAQSVGAVL VATWIHDYAP AIRGLVLASP AFKVKLYVPL ARPALALWHR LRGLFFINSY VKGRYLTHDR QRGASFNDP LITRAIAVNI LLDLYKTSER IIRDAAAITL PTQLLISGDD YVHRQPQID FYQRLRSPLK ELHLLPGFYH DTLGEENRAL AFEKMQSFIS RLYANKSQKF DYQHEDCTGP SADRWRLLSG GPVPLSPVDL AYRFMRKAMK LFGTHSSGLH LGMSTGFDSG SSLDYVYQNG PQGSNAFGRL VDKIYLN SVG WRGIRQRKTH LQILIKQAVA DLHAKGLAVR VVDIAAGHGR YVLDALANEP AVSDILLRDY SELNVAQQQE MIAQRGMSGR VRFEQGDAFN PEELSALTPR PTLAIVSGLY ELFPENEQVK NSLAGLANAI EPGGILIYTG QPWHPQLEMI AGVLTSHKDG KPWVMRVR SQ GEMDSLVRDA GFDKCTQRID EWGIFTVSMA VRRDN
Source	E.coli
Target Names	ynbC
Protein Names	Recommended name: Uncharacterized protein ynbC
Expression Region	1-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.