



# Recombinant Escherichia coli Uncharacterized protein ynbC (ynbC)

<b>Product Code</b>	CSB-MP300535ENV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P76092
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MENSRIPGEH FFTTSDNTAL FYRHWPALQP GAKKVIVLFH RGHEHSGRLQ HLVDELAMPD TAFYAWDARG HGKSSGPRGY SPLSARSVRD VDEFVRF AAS DSQVGL EEVV VIAQSVGAVL VATWIHDYAP AIRGLVLASP AFKVKLYVPL ARPALALWHR LRGLFFINSY VKGRYLTHDR QRGASFNDP LITRAIAVNI LLDLYKT SER IIRDAAAITL PTQL LISGDD YVHRQPQID FYQRLRSPLK ELHLLPGFYH DTLGEENRAL AFEKMQSFIS RLYANKSQKF DYQHEDCTGP SADRWRLLSG GPVPLSPVDL AYRFMRKAMK LFGTHSSGLH LGMSTGFDSG SSLDYVYQNQ PQGSNAFGRL VDKIYLN SVG WRGIRQRKTH LQILIKQAVA DLHAKGLAVR VVDIAAGHGR YVLDALANEP AVSDILLRDY SELNVAQQQE MIAQRGMSGR VRFEQGD AFN PEELSALTPR PTLAIVSGLY ELFPENEQVK NSLAGLANAI EPGGILIYTG QPWHPQLEMI AGVLTSHKDG KPWVMRVR SQ GEMDSLVRDA GFDKCTQRID EWGIFTVSMA VRRDN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ynbC
<b>Protein Names</b>	Recommended name: Uncharacterized protein ynbC
<b>Expression Region</b>	1-585
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.