



Recombinant Escherichia coli Probable diguanylate cyclase YeaP (yeaP)

Product Code	CSB-BP300549ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76245
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSDQIIARVS QSLAKEQSLE SLVRQLLEML EMVTDMESTY LTKVDVEARL QHIMFARNSQ KMYIPENFTV SWDYSLCKRA IDENCFFSDE VPDRWGDICIA ARNLGITTFLL STPIHLPDGS FYGTLCAASS EKRQWSERAQ QVLQLFAGLI AQYIQKEALV EQLREANAAL IAQSYTDSLTL GLPNRRAIFE NLTTLFSLAR HLNHKIMIAF IDLDNFKLIN DRFGHNSGDL FLIQVGERLN TLQQNGEVIG RLGGEFLVV SLNNENADIS SLRERIQQQI RGEYHLGDVD LYYPGASLGI VEVDPETTDA DSALHAADIA MYQEKKHKQK TPFVAHPALH S
Source	Baculovirus
Target Names	dgcP
Protein Names	Recommended name: Probable diguanylate cyclase YeaP Short name= DGC EC= 2.7.7.65
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.