



Recombinant Escherichia coli Autoinducer 2 import ATP-binding protein LsrA (IsrA)

Product Code	CSB-BP300580ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77257
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MQTS DTRALP LLCARSVYKQ YSGVNV LKGI DFTLHQGEVH ALLGGNGAGK STLMKIIAGI TPADSGTLEI EGN NYVRLTP VHAHQ LGIYL VPQEPLLFPS LSIKENILFG LAKKQLSMQK MKNLLAALGC QFDLHSLAGS LDVADRQMVE ILRGLMRDSR ILILDEPTAS LTPAETERLF SRLQELLATG VGIVFISHKL PEIRQIADRI SVMRDGTIAL SGKTSELSTD DIIQAITPAV REKSL SASQK LWLELPGNRP QHAAGTPVLT LENLTGEGFR NVSLTLNAGE ILGLAGLVGA GRTELAETLY GLRTL RGGRI MLNGKEINKL STGERLLRGL VYLPEDRQSS GLNLDASLAW NVCALTHNLR GFWAKTAKDN ATLERYRRAL NIKFNQPEQA ARTLSGGNQQ KILIAKCLEA SPQVLIVDEP TRGVDVSARN DIYQLLR SIA AQNVA VLLIS SDLEEIELMA DRVYVMHQGE ITHSALTERD INVETIMRVA FGDSQRQEAS C</p>
Source	Baculovirus
Target Names	IsrA
Protein Names	Recommended name: Autoinducer 2 import ATP-binding protein LsrA Short name= AI-2 import ATP-binding protein LsrA EC= 3.6.3.- Alternative name(s): EGO10A
Expression Region	1-511
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.