



Recombinant Escherichia coli Uncharacterized J domain-containing protein djIC (djIC)

Product Code	CSB-BP300586ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77359
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKTCWQILEI ESTTQIDIIR QAYLARLPLC HPETDPQGFK ALRQAYEEAL RLAVNPVEEA DDEEKDAAAE HEILRAFRTL LDESDFRFP SAWQKFIQQL NTWNMEDVDQ LRWPLCAIAI EARYLSLNCA SLLAERLNWH SFNDSEGMDE EEREAFLEAI QAGDCDFLS LLEYPIALQN QTVEYYFALE RCCRYHPDYV TAFLAMEGPW LIPDDAKLHR KLLRWYSSVQ TGMAELIPVA QQWQTEEPES EDARYYLCAQ RLYCGEGESL LADLCAYWES YPSTQADNLL LQWSKRHCPD YFALLVMVIE ARSMVDAQGQ PLKYVPGESA RTRLLWAEIL HSGKLSPLGQ SFIESLFFKR KAWAWWKS RV GSETEQDSPF LDLYRVAEQV VLEAFPQEM LARLNRLEG GDAHPLEAIV TRMLLTKVKL EPEDEDVDEP TPENHEEKND EGEKPQSITS IIKISLTVLV IGYALGKIAM LFS
Source	Baculovirus
Target Names	djIC
Protein Names	Recommended name: Uncharacterized J domain-containing protein djIC
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.