



Recombinant Escherichia coli Uncharacterized fimbrial-like protein ydeQ (ydeQ)

Product Code	CSB-MP300599ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77588
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	CNVD GGSSIGAGTT SVYVNLDPVI QPGQNLVVDL SQHISCWNDY GGWYDTDHIN LVQGSFAFAGS LQSYKGSLYW NNVTYPFPLT TNTNVLDIGD KTPMPLPLKL YITPVGAAGG VVIKAGEVIA RIHMYKIATL GSGNPRNFTW NIISNNVVM PTGGCTVDSR NVTVDLPDFP GSAEIPLGVY CSSEQKLSFY LSGATTDSSR QVFANTAPDA TKASGVGVTL MRNGKILATG ENVSLGTVNK SKVPLGLSAT YGQTGNKVSA GTVQSVIGVT FIYE
Source	Mammalian cell
Target Names	ydeQ
Protein Names	Recommended name: Uncharacterized fimbrial-like protein ydeQ
Expression Region	27-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.