



Recombinant Escherichia coli FeS cluster assembly protein sufD (sufD)

Product Code	CSB-BP300604ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77689
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MAGLPNSSNA LQQWHHLFEA EGTKRSPQAAQ QHLQQLLRTG LPTRKHENWK YTPLEGLINS QFVSIAGEIS PQRDALALT LDSVRLVFVD GRYVPALSDA TEGSGYEVSI NDDRQGLPDA IQAEVFLHLT ESLAQSVTHI AVKRGQRPAAK PLLLMHITQG VAGEEVNTAH YRHHLDLAEG AEATVIEHFV SLNDARHFTG ARFTINVAAN AHLQHIKLAFL ENPLSHHFAH NDLLAEDAT AFSHSFLGG AVLRHNTSTQ LNGENSTLRI NSLAMPVKNE VCDTRTWLEH NKGFCNSRQL HKTIVSDKGR AVFNGLINVA QHAIKTDGQM TNNLLMGKL AEVDTKPQLE IYADDVKCSH GATVGRIDDE QIFYLRSRGI NQQDAQQMII YAFAAELTEA LRDEGLKQQV LARIGQRLPG GAR
Source	Baculovirus
Target Names	sufD
Protein Names	Recommended name: FeS cluster assembly protein sufD
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.