



Recombinant Escherichia coli Propionate catabolism operon regulatory protein (prpR)

Product Code	CSB-BP300607ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77743
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MAHPPRLNDD KPVIWTVSVT RLFELFRDIS LEFDHLANIT PIQLGFEEKAV TYIRKKLANE RCDAlIAAGS NGAYLKSRLS VPVILIKPSG YDVLQALAKA GKLTSSIGVV TYQETIPALV AFQKTFNLRL DQRSYITEED ARGQINELKA NGTEAVVGAG LITDLAEEAG MTGIFIYSAA TVRQAFSDAL DMTRMSLRHN THDATRNALR TRYVLGDMLG QSPQMEQVRQ TILLYARSSA AVLIEGETGT GKELAAQAIH REYFARHDAR QGKKSHPFVA VNCGAI AESL LEAELFGYEE GAFTGSRGG RAGLFEIAHG GTLFLDEIGE MPLPLQTRLL RVLEEKEVTR VGGHQVPVD VRVISATHCN LEEDMQQGRF RRDLFYRLSI LRLQLPPLRE RVADILPLAE SFLKVSLAAL SAPFSAALRQ GLQASETVLL HYDWPGNIRE LRNMMERLAL FLSVEPTPDL TPQFMQLLLP ELARES AKTP APRLLTPQQA LEKFNGDKTA AANYLGISRT TFWRRLKS
Source	Baculovirus
Target Names	prpR
Protein Names	Recommended name: Propionate catabolism operon regulatory protein
Expression Region	1-528
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.