



Recombinant Escherichia coli Uncharacterized HTH-type transcriptional regulator yhaJ (yhaJ)

Product Code	CSB-EP300687ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67660
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MAKERALTLE ALRVMDAIDR RGSFAAAADE LGRVPSALSY TMQKLEEELD VVLFD RSGHR TKFTNVGRML LERGRVLEA ADKLT TDAEA LARGWETHLT IVTEALVPTP AFFPLIDKLA AKANTQLAII TEVLAGAWER LEQGRADIVI APDMHFRSSS EINSRKLYTL MNVYVAAPDH PIHQEPEPLS EVTRVKYRGI AVADTARERP VLT VQLLDKQ PRLTVSTIED KRQALLAGLG VATMPYPMVE KDIAEGRLRV VSPESTSEID IIMAWRRDSM GEAKSWCLRE IPKLFNGK
Source	E.coli
Target Names	yhaJ
Protein Names	Recommended name: Uncharacterized HTH-type transcriptional regulator yhaJ
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.