



Recombinant Rauvolfia serpentina Strictosidine synthase (STR1)

Product Code	CSB-BP300722RBY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68175
Product Type	Recombinant Protein
Immunogen Species	Rauvolfia serpentina (Serpentine wood) (Ophioxylon serpentinum)
Purity	≥85% (SDS-PAGE)
Sequence	SLALSSPI LKEILIEAPS YAPNSFTFDS TNKGFYTSVQ DGRVIKYEGP NSGFVDFAYA SPYWNKAFCE NSTDAEKRPL CGRTYDISYN LQNNQLYIVD CYYHLSVVGS EGGHATQLAT SVDGVPFKWL YAVTVDQRTG IVYFTDVSTL YDDRGVQQIM DTSDKTGRLI KYDPSTKETT LLLKELHVP GAEVSADSSF VLVAEFLSHQ IVKYWLEGPK KGTAEVLVKI PNPNGNIKRNA DGHFVWSSSE ELDGNMHGRV DPKGKIFDEF GNILEVIPLP PPFAGEHFEQ IQEHDGLLYI GTLFHGSGVI LVIYDKKGN SF VSSH
Source	Baculovirus
Target Names	STR1
Protein Names	Recommended name: Strictosidine synthase EC= 4.3.3.2
Expression Region	23-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.