



Recombinant Arabidopsis thaliana Dof zinc finger protein DOF1.2 (DOF1.2)

Product Code	CSB-MP300730DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68349
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MLPYIGHNSY QQHQFPLPEM EIPEKWKLSY EQEAITAPAC PRCASSNTKF CYYNNYSLSQ PRYFCKGCRR YWTKGGSLRN IPVGGGCRKR SRSRQNSHKR FGRNENRPDG LINQDDGFQS SPPGSDIDLA AVFAQYVTDR SPSSTDNTTG SDQDSPITTT THALESLSWD ICQETDVDLG FYGEFNNLTQ KTKEDQEVFG QFLQEDREEI FEFQGLLDDK EIQEILECSF SEEPDQLVSQ GSFMINGDNW SSTDLTRFGI
Source	Mammalian cell
Target Names	DOF1.2
Protein Names	Recommended name: Dof zinc finger protein DOF1.2 Short name= AtDOF1.2
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.