



Recombinant Bacillus subtilis Phosphate acyltransferase (plsX)

Product Code	CSB-EP300932BRJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P71018
Product Type	Recombinant Protein
Immunogen Species	Bacillus subtilis (strain 168)
Purity	>85% (SDS-PAGE)
Sequence	MRIAVDAMGG DHAPKAVIDG VIKGIEAFDD LHITLVGDKT TIESHLTTTS DRITVLHADE VIEPTDEPVR AVRRKKNSSM VLMAQEVAEN RADACISAGN TGALMTAGLF IVGRIKGIDR PALAPTLPTV SGDGFLLLDV GANVDAKPEH LVQYAIMGSV YSQQVRGVTs PRVGLLNVGT EDKKGNELTK QTFQILKETA NINFIGNVEA RDLLDDVADV VVDGFTGNV TLKTLEGSAL SIFKMMRDVM TSTLTSKLAA AVLKPKLKEM KMKMEYSNYG GASLFGLKAP VIKAHGSSDS NAVFHAIRQA REMVSQNVAA LIQEEVKEEK TDE
Source	E.coli
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]-phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.