



Recombinant *Synechococcus elongatus* tRNA (guanine-N (7)-)-methyltransferase

Product Code	CSB-BP300979FPY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P72546
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	≥85% (SDS-PAGE)
Sequence	MARVRVRQHV NPLSQKFQVV TTWPDWQQVY ADCDRPLHLD IGCARGRFL AMATAQPEWN YLGLIUREPL VDEANAIARE RELTNLYYHF SNANLDLEPL LRS�PTGILQ RVSİQFPDPW FKRRHQKRRV VQPELVQALA TALPAGAEVF LQSDVLEVQA EMCEHFAAEP RFQRTCLDWL PENPLPVPT REIAVQNKQL PVYRALFIRQ PAD
Source	Baculovirus
Target Names	trmB
Protein Names	Recommended name: tRNA (guanine-N(7)-)-methyltransferase EC= 2.1.1.33 Alternative name(s): tRNA (guanine(46)-N(7))-methyltransferase tRNA(m7G46)- methyltransferase
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.