



# Recombinant Macaca mulatta (Rhesus macaque) Synaptosomal-associated protein 25 (SNAP25)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-MP021873   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.  |
| <b>Uniprot No.</b>       | P60877   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Macaca mulatta (Rhesus macaque)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MAEDADMRNELEEMQRRADQLADESLESTRMLQLVEESKDAGIRTLVMLDE<br>QGEQLERI<br>EEGMDQINKDMKEAEKNLTDLGKFCGLVCPCNKLKSSDAYKKAWGNNQDG<br>VVASQPARV<br>VDEREQMAISGGFIRRVVTNDARENEMDENLEQVSGIIGNLRHMALDMGNEIDT<br>QNRQIDR IMEKADSNKTRIDEANQRATKMLGSG   |
| <b>Source</b>            | Mammalian cell   |
| <b>Target Names</b>      | SNAP25   |
| <b>Protein Names</b>     | Recommended name: Synaptosomal-associated protein 25 Short name= SNAP-25 Alternative name(s): Synaptosomal-associated 25 kDa protein   |
| <b>Expression Region</b> | 1-206  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full length protein  |
| <b>Target Details</b>    | Synaptic vesicle membrane docking and fusion is mediated by SNAREs (soluble N-ethylmaleimide-sensitive factor attachment protein receptors) located on the vesicle membrane (v-SNAREs) and the target membrane (t-SNAREs). The assembled v-SNARE/t-SNARE complex consists of a bundle of four helices, one of which is supplied by v-SNARE and the other three by t-SNARE. For t-SNAREs on the plasma membrane, the protein syntaxin supplies one helix and This protein contributes the other two. Therefore, this gene product is a presynaptic plasma membrane protein involved in the regulation of neurotransmitter release. Two alternative transcript variants encoding different protein isoforms have been described for this gene. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.