



Recombinant Mouse NEDD8-conjugating enzyme Ubc12 (Ube2m)

Product Code	CSB-BP025468MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61082
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MIKLFSLKQQ KKEEESAGGT KGSSKKASAA QLRIQKDINE LNLPKTCDIS FSDPDDLINF KLVICPDEGF YKSGKFVFSF KVGQGYPHDP PKVKCETMVY HPNIDLEGNV CLNILREDWK PVLINSIIY GLQYLFLEPN PEDPLNKEAA EVLQNNRRLF EQNVQRSMRG GYIGSTYFER CLK
Source	Baculovirus
Target Names	Ube2m
Protein Names	Recommended name: NEDD8-conjugating enzyme Ubc12 EC= 6.3.2.- Alternative name(s): NEDD8 carrier protein NEDD8 protein ligase Ubiquitin-conjugating enzyme E2 M
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein is linked with a ubiquitin-like protein, NEDD8, which can be conjugated to cellular proteins, such as Cdc53/culin.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.