



Recombinant Escherichia coli 2-octaprenyl-3-methyl-6-methoxy-1,4-benzoquinol hydroxylase (ubiF)

Product Code	CSB-MP301152ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75728
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTNQPTETIAI VGGGMVGGAL ALGLAQHGFA VTVIEHAEPA PFVADSQPDV RISAISAASV SLLKGLGVWD AVQAMRCHPY RRLETWEWET AHVVFDAEEL KLPLLLGYMVE NTVLQQALWQ ALEAHPKVTL RVPGLIALH RHDDLQELEL KGGEVIRAKL VIGADGANSQ VRQMAGIGVH AWQYAQSCML ISVQCENDPG DSTWQQFTPD GPRAFLPLFD NwasLVWYDS PARIRQLQNM NMAQLQAEIA KHFPSRLGYV TPLAAGAFPL TRRHALQYVQ PGLALVGDAA HTIHPLAGQG VNLGYRDVDA LIDVLVNARS YGEAWASYPV LKRYQMRRMA DNFIMQSGMD LFYAGFSNNL PPLRFMRNLG LMAAERAGVL KRQALKYALG L
Source	Mammalian cell
Target Names	ubiF
Protein Names	Recommended name: 2-octaprenyl-3-methyl-6-methoxy-1,4-benzoquinol hydroxylase EC= 1.14.13.-
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.