



Recombinant Escherichia coli Uncharacterized protein ybjT (ybjT)

Product Code	CSB-BP301158ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75822
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MPQRILVLGA SGYIGQHLVR TLSQQGHQIL AAARHVDRLA KLQLANVSCH KVDLSWPDNL PALLQDIDTV YFLVHSMGEG GDFIAQERQV ALNVRDALRE VPVKQLIFLS SLQAPPHEQS DHLRARQATA DILREANVPV TELRAGIIVG AGSAAFEVMR DMVYNLPVLT PPRWVRSRTT PIALENLLHY LVALLDHPAS EHRIFEAAGP EVLSYQQQFE HFMAVSGKRR WLIPIPLPTR WISVWFLNVI TSVPPTTARA LIQGLKHDLL ADDTALRALI PQRLIAFDDA VRSTLKEEEK LVNSSDWGYD AQAFARWRPE YGYFAKQAGF TVKTSASLAA LWQVVNQIGG KERYFFGNIL WQTRALMDRA IGHKLAKGRP EREYLQTGDA VDSWKVIVVE PEKQLTLLFG MKAPGLGRLC FSLEDKGDYR TIDVRAFWHP HGMPGLFYWL LMIPAHLFIF RGMAKQIARL AEQSTD
Source	Baculovirus
Target Names	ybjT
Protein Names	Recommended name: Uncharacterized protein ybjT
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.