



# Recombinant Escherichia coli Uncharacterized protein ybjT (ybjT)

<b>Product Code</b>	CSB-YP301158ENV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P75822
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPQRILVLGA SGYIGQHLVR TLSQQGHQIL AAARHVDRLA KLQLANVSCH KVDLSWPDNL PALLQDIDTV YFLVHSMGEG GDFIAQERQV ALNVRDALRE VPVKQLIFLS SLQAPPHEQS DHLRARQATA DILREANVPV TELRAGIIVG AGSAAFEVMR DMVYNLPVLT PPRWVRSRTT PIALENLLHY LVALLDHPAS EHRIFEAAGP EVLSYQQQFE HFMAVSGKRR WLIPIPLPTR WISVWFLNVI TSVPPTTARA LIQGLKHDLL ADDTALRALI PQRLIAFDDA VRSTLKEEEK LVNSSDWGYD AQAFARWRPE YGYFAKQAGF TVKTSASLAA LWQVVNQIGG KERYFFGNIL WQTRALMDRA IGHKLAKGRP EREYLQTGDA VDSWKVIVVE PEKQLTLLFG MKAPGLGRLC FSLEDKGDYR TIDVRAFWHP HGMPGLFYWL LMIPAHLFIF RGMAKQIARL AEQSTD
<b>Source</b>	Yeast
<b>Target Names</b>	ybjT
<b>Protein Names</b>	Recommended name: Uncharacterized protein ybjT
<b>Expression Region</b>	1-476
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.