



Recombinant Human Actin-related protein 2 (ACTR2)

Product Code	CSB-MP001248HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61160
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSQGRKVVV CDNGTGFVKC GYAGSNFPEH IFPALVGRPI IRSTTKVGNIEIKDLMVGDE ASELRSMLEV NYPMENGIVR NWDDMKHLWD YTFGPEKLNIDTRNCKILLT EPPMNPTKNR EKIVEVMFET YQFSGVYVAI QAVLTLYAQQLLTGVVVDSG DGVTHICPVY EGFSLPHLTR RLDIAGRDI RYLIKLLLLR GYAFNHSADF ETVRMIKEKL CYVGYNIEQE QKLALETTVL VESYTLDPGR IIKVGGERFE APEALFQPHL INVEGVGVAE LLFNTIQAAD IDTRSEFYKH IVLSGGSTMY PGLPSRLERE LKQLYLERVL KGDVEKLSKF KIRIEDPPRR KHMVFLGGAV LADIMKDKDN FWMTRQEYQE KGVRVLEKLG VTVR
Source	Mammalian cell
Target Names	ACTR2
Protein Names	Recommended name: Actin-related protein 2 Alternative name(s): Actin-like protein 2
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.