



Recombinant Escherichia coli Uncharacterized protein ycbX (ycbX)

Product Code	CSB-EP301162ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75863
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MATLIRLFIH PVKSMRGIGL THALADVSGL AFDRIFMITE PDGTFITARQ FPQMVRFTPS PVHDGLHLTA PDGSSAYVRF ADFATQDAPT EVWGTHFTAR IAPDAINKWL SGFFSREVQL RWVGPQMTRR VKRHNTVPLS FADGYPYLLA NEASLRDLQQ RCPASVKMEQ FRPNLVVSGA SAWEEDRWKV IRIGDVVFDV VKPCSRCIFT TVSPEKGQKH PAGEPLKTLQ SFRTAQDNGD VDFGQNLIAR NSGVIRVGDE VEILATAPAK IYGAAAADDT ANITQQPDAN VIDIDWQQQAF RGNNQQVLE QLENQGIRIP YSCRAGICGS CRVQLLEGEV TPLKKSAMGD DGTILCCSCV PKTALKLAR
Source	E.coli
Target Names	ycbX
Protein Names	Recommended name: Uncharacterized protein ycbX
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.