



Recombinant Escherichia coli Uncharacterized protein yfaQ (yfaQ)

Product Code	CSB-MP301203ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76463
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	E ETPLQLVLRG AQHDQLYQLS SSGVTKVSAL PDSLTTPLGS LWKLYVYAWL EDTHQPEQPY QCRGNSPEEV YCCQAGESIT RDTALVRSCG LYFAPQRLHI GADVWGQYWQ QRQAPAWLAS LTTLKPETSV TVKSLDLSLA TLPAQNKAQE VLLDVVLDEA KIGVASMLGS RVRVKTWSWF ADDKQEIRQG GFAGWLT DGT PLWVTGSGTS KTVLTRYATV LNRVLPVPTQ VASGQCVEVE LFARYPLKKI TAEKSTTAVN PGVLNGRYRV TFTNGNHITF VSHGETTLLS EKGKLLQSH LDREEYVARV LDREAKSTPP EAAKAMTVAI RTFLQQNANR EGDCLTIPDS SATQRVSASP ATTGARTMTA WTQDLIYAGD PVHYHGS RAT EGTLSWRQAT AQAGQGERYD QILAFAYPDN SLSRWGAPRS TCQLLPKAKA WLAKKMPQWR RILQAETGYN EPDVFAVCRL VSGFPYTDRQ QKRLFIRNFF TLQDRLDLTH EYLHLAFDGY PTGLDENYIE TLTRQLLMD
Source	Mammalian cell
Target Names	yfaQ
Protein Names	Recommended name: Uncharacterized protein yfaQ
Expression Region	20-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.