



Recombinant Escherichia coli Aminopeptidase ypdF (ypdF)

Product Code	CSB-YP301207ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76524
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MTLLASLRDW LKAQQQLDAVL LSSRQNKQPH LGISTGSGYV VISRESAHIL VDSRYYYEVE ARAQGYQLHL LDATNTLTTI VNQIIADEQL QTLGFEGQQV SWETAHRWQS ELNAKLVSAT PDVLRQIKTP EEVEKIRLAC GIADRGAEHI RRFIQAGMSE REIAAELEWF MRQQGAEKAS FDTIVASGWR GALPHGKASD KIVAAGEFVT LDFGALYQGY CSDMTRTLLV NGEVSAESH LLFNVYQIVL QAQLAAISAI RPGVRCQQVD DAARRVITEA GYGDYFGHNT GHAIGIEVHE DPRFSPRDTT TLQPGMLLTV EPGIYLPQQG GVRIEDVVLV TPQGAEVLYA MPKTVLLTGE A
Source	Yeast
Target Names	ypdF
Protein Names	Recommended name: Aminopeptidase ypdF EC= 3.4.11.-
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.