



Recombinant Escherichia coli FeS cluster assembly protein sufB (sufB)

Product Code	CSB-BP301237ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77522
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSRNTEATDD VKTWTGGPLN YKEGFFTQLA TDELAKGINE EVVRAISAKR NEPEWMLEFR LNAYRAWLEM EEPHWLKAHY DKLNYQDYSY YSAPSCGNCD DTCASEPGAV QQTGANAFSL KEVEAAFEQL GVPVREGKEV AVDAIFDSVS VATTYREKLA EQGIIFCSFG EAIHDHPELV RKYLGTVVPG NDNFFAALNA AVASDGTFIY VPKGVRCPME LSTYFRINAE KTGQFERTIL VADEDSYVSY IEGCSAPVRD SYQLHAAVVE VIIHKNAEVK YSTVQNWFPG DNNTGGILNF VTKRALCEGE NSKMSWTQSE TGSAITWKYP SCILRGDNSI GEFYSVALTS GHQQADTGTK MIHIGKNTKS TIISKGISAG HSQNSYRGLV KIMPTATNAR NFTQCDSMLI GANCGAHTFP YVECRNNSAQ LEHEATTSRI GEDQLFYCLQ RGISEEDAIS MIVNGFCKDV FSELPLEFAV EAQKLLAISL EHSV G
Source	Baculovirus
Target Names	sufB
Protein Names	Recommended name: FeS cluster assembly protein sufB
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.