



Recombinant Schizosaccharomyces pombe Putative ribosomal RNA methyltransferase C2G2.15c (SPBC2G2.15c)

Product Code	CSB-BP301276SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P78860
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MFGSKTLFSW AAKRSKDFYR KSKIDNFRS RAAYKLIELN SKYRFINKED VVIDVGFAPG SWSQVAKKLV GNKGKVIDIG IQHIAPPEGV LPIYGDIRD NTLTKLFEAL RLLHEPNTND SIDCRVVDVAV ISDMLHKATG IRIRDHALSM ELCASALHVA LTFLKSNNGSF ICKFYMGDED ADLQNLKSH FRFVQVMKPK ASLKESREAY FVCLERKP
Source	Baculovirus
Target Names	mrm2
Protein Names	Recommended name: Putative ribosomal RNA methyltransferase C2G2.15c EC= 2.1.1.- Alternative name(s): rRNA (uridine-2'-O-)-methyltransferase
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.