



# Recombinant Physcomitrella patens subsp. patens Myb-related protein Pp1 (PP1)

<b>Product Code</b>	CSB-EP301296EVD
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P80074
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Physcomitrella patens subsp. patens (Moss)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LGNRWSAIAI PRRTDNEIKN YWNTHLKKRL MQMGIDPVTH KSTAAEELVH YSIIPGLRPV VSTNLTHMSQ WDSARAEAEA RLSRQSSLTS PASDLAQTSL EHQKSNLSTK NDVNNQVASS NFMSSWKAQV TETLRPNFGV VELDKPPASP VNLQKFLQEW ESSLKAPQPE MEGPSIHDSI TNIPSLSSGT ASELVSTQYS PDVSSASYRM PVSSMILPSP AQISERLMAS RHCDSLLPLP RIDLAGSNIF PAQSPSATES SFSFSGGLCKV DNEHQYSPTS ILHGPSGHDS SYSSPCGSSS SAYDSVDILA QTFSLVHNDN HQHSRATLAM NFSQEPSFWT QQQVALVDAL QPESYFTHL GVNIAQKPPS HSFNLSVNSM VNQVGVPNNI PFHSLICVE SIFVAYFRRS RCTSIANLHM FFVFRFYFYL FYFSFDFYRV PLRRRKSQW
<b>Source</b>	E.coli
<b>Target Names</b>	PP1
<b>Protein Names</b>	Recommended name: Myb-related protein Pp1
<b>Expression Region</b>	1-449
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.