



Recombinant *Thermus thermophilus* 30S ribosomal protein S2 (rpsB)

Product Code	CSB-BP301325TNT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P80371
Product Type	Recombinant Protein
Immunogen Species	<i>Thermus thermophilus</i> (strain HB8 / ATCC 27634 / DSM 579)
Purity	>85% (SDS-PAGE)
Sequence	PVEITVKEL LEAGVHFGHE RKRWNPKFAR YIYAERNGIH IIDLQKTMEELERTFRFIED LAMRGGTILF VGTKKQAQDI VRMEAERAGM PYVNQRWLGGMLTNFKTISQ RVHRLEELEA LFASPEIEER PKKEQVRLKH ELERLQKYLSGFRLLKRLPD AIFVVDPTKE AIAVREARKL FIPVIALADT DSDPDLVDYIIPGNDDAIRS IQLILSRAVD LIIQARGGVV EPSPSYALVQ EAEATETPEG ESEVEA
Source	Baculovirus
Target Names	rpsB
Protein Names	Recommended name: 30S ribosomal protein S2
Expression Region	2-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.