



Recombinant Escherichia coli Alkanesulfonate monooxygenase (ssuD)

Product Code	CSB-EP301373ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P80645
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	SLNMFWFLP THGDGHYLG T EEGSRPVDHG YLQQIAQAAD RLG YTGVLIP TGRSCEDAWL VAASMIPVTQ RLKFLVALRP SVTSPTVAAR QAATLDRLSN GRALFNLVTG SDPQELAGDG VFLDHSERYE ASAEFTQVWR RLLQRETVD F NGKHIHVRGA KLLFPAIQQP YPLYFGGSS DVAQELAAEQ VDLYLTWGE P PELVKEKIEQ VRAKAAAHGR KIRFGIRLHV IVRETNDEAW QAAERLISHL DDETIAKAQA AFARTDSVGQ QRMAALHNGK RDNLEISPNL WAGVGLVRGG AGTALVGDGP TVAARINEYA ALGIDSFVLS GYPHLEEAYR VGELLFPLLD VAIPEIQPQ PLNPQGEAVA NDFIPRKVAQ S
Source	E.coli
Target Names	ssuD
Protein Names	Recommended name: Alkanesulfonate monooxygenase EC= 1.14.14.5 Alternative name(s): FMNH ₂ -dependent aliphatic sulfonate monooxygenase Sulfate starvation-induced protein 6 Short name= SSI6
Expression Region	2-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.