



Recombinant Human WD repeat-containing protein 5 (WDR5)

Product Code	CSB-BP026038HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61964
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ATEEKKPET EAARAQPTPS SSATQSKPTP VKPNYALKFT LAGHTKAVSS VKFSPNGEWL ASSADKLIK IWGAYDGKFE KTISGHKLG I SDVAWSSDSN LLVSASDDKT LKIWDVSSGK CLKTLKGHSN YVFCCNFNPQ SNLIVSGSFD ESVRIWDVKT GKCLKTLP AH SDPVSAVHFN RDGSLIVSSS YDGLCRIWDT ASGQCLKTLI DDDNPPVSFV KFSPNGKYIL AATLDNTLKL WDYSKGKCLK TYTGHKNEKY CIFANFSVTG GKWIVSGSED NLVYIWNLQT KEIVQKLQGH TDVVISTACH PTENIIASAA LENDKTIKLW KSDC
Source	Baculovirus
Target Names	WDR5
Protein Names	Recommended name: WD repeat-containing protein 5 Alternative name(s): BMP2-induced 3-kb gene protein
Expression Region	2-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation. This protein contains 7 WD repeats. Alternatively spliced transcript variants encoding the same protein have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.