



Recombinant Enterobacteria phage lambda Protein rexA (rexA)

Product Code	CSB-BP301608ECW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68924
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Purity	>85% (SDS-PAGE)
Sequence	MKNGFYATYR SKNKGKDKRS INLSVFLNSL LADNHHLQVG SNLYIHKID GKTFLFTKTN DKSLVQKINR SKASVEDIKN SLADDES LGF PSFLFVEGDT IGFARTVFGP TTSDLTDFLI GKGMSLSSGE RVQIEPLMRG TTKDDVMHMH FIGRTTVKVE AKLPVFGDIL KVLGATDIEG ELFDSLDIVI KPKFKRDIKK VAKDIIFNPS PQFSDISLRA KDEAGDILTE HYLSEKGHLS APLNKVTNAE IAEEMAYCYA RMKSDILECF KRQVGKVKD
Source	Baculovirus
Target Names	rexA
Protein Names	Recommended name: Protein rexA
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.