



Recombinant *Geobacillus stearothermophilus* Type-2 restriction enzyme BsoBI (bsoBIR)

Product Code	CSB-BP301752GFM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70985
Product Type	Recombinant Protein
Immunogen Species	<i>Geobacillus stearothermophilus</i> (<i>Bacillus stearothermophilus</i>)
Purity	>85% (SDS-PAGE)
Sequence	MNTQKPFENH LKSVDDLKTT YEEYRAGFIA FALEKNKRST PYIERARALK VAASVAKTPK DLLYLEDIQD ALLYASGISD KAKKFLTEDD KKESINNLIE NFLEPAGEEF IDELIFRYLL FQGDSLGGTM RNIAGALAQQ KLTRAIISAL DIANIPYKWL DSRDKKYTNW MDKPEDDYEL ETFAKGISWT INKGHRTLMY NITVPLVKKV VDICLFNCEP EIYTPQKVHQ QPEKYLLLGE LKGGIDPAGA DEHWKTANTA LTRIRNKFSE KGLSPKTIFI GAAIEHSMAE EIWDQLQSGS LTNSANLTKT EQVGS�CRWI INI
Source	Baculovirus
Target Names	bsoBIR
Protein Names	Recommended name: Type-2 restriction enzyme BsoBI Short name= R.BsoBI EC= 3.1.21.4 Alternative name(s): Endonuclease BsoBI Type II restriction enzyme BsoBI
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.