



# Recombinant Haemophilus influenzae Trk system potassium uptake protein trkA (trkA)

<b>Product Code</b>	CSB-MP301768HTA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P71354
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Haemophilus influenzae (strain ATCC 51907 / DSM 11121 / KW20 / Rd)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKIIILGAGQ VGTTLAENLV SEDNDITLVD NESPQLQTLQ EKHDLRVVQG SPSSPKVLRD AGAADADLMV AVTASDEINM VACQMGYTLF NTPTRIARIR NSEYLREKDK LFNNENIPID HLISPENLVT DEITRLIAYP GALQVAHFAN NRISIVVKA YYGGALVGYA LSAFREHMPH IDCRIMSILR NGKPIRPQGS TIVEAGDEIT FICATEHIKA IMGELQRLEK PYKRVMIVGG GNAVFGVAKR LENSCTVKLI ERDSNRAQAL AEKLPKTLVF NGDASDQNLL FEEHIESVDV FLSLSSDDEA NIMSALLAKR LGAKKAMVLI QRIAYINLIQ GGTIDIAVSP QQVTISALLG HVRKGDVKNV ATRLRHGIAEA IEIVAHGNVN TSNIVGRKIG ELRLPMGIII GALLRGNDVI IARRQVIIIE GDHIVIYLSLSD KKNVPEIEKL FQPSAFFI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	trkA
<b>Protein Names</b>	Recommended name: Trk system potassium uptake protein trkA Short name=K(+)-uptake protein trkA
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.