



Recombinant Human 40S ribosomal protein S6 (RPS6)

Product Code	CSB-BP020463HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P62753
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKLNISFPAT GCQKLIEVDD ERKLRIFYEK RMATEVAADA LGEEWKGYYV RISGGNDKQG FPMKQGVLTG GRVRLLSKSG HSCYRPRRTG ERKRKSVRGC IVDANLSVLN LVIVKKGEKD IPGLTDTTVP RRLGPKRASR IRKLFNLSKE DDVRQYVVRK PLNKEGKKPR TKAPKIQRLV TPRVLQHKRR RIALKKQRTK KNKEEAAEYA KLLAKRMKEA KEKRQEQIAK RRRLSSLRAS TSKSESSQK
Source	Baculovirus
Target Names	RPS6
Protein Names	Recommended name: 40S ribosomal protein S6 Alternative name(s): Phosphoprotein NP33
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.