



# Recombinant Mycoplasma pneumoniae Putative type-1 restriction enzyme specificity protein MPN\_290 (MPN\_290)

<b>Product Code</b>	CSB-EP301952MLW-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P75487
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVYSKTFRVE EKSITVSARG TIGVVFYRDF AYLPAVSLIC FVPKEEFDIR FLFHALRAIK FKKQGSATGQ LTVAQFKEYG IHVPSLKKKQK EIAAILDPLY SFFTDLNEGL PAEIELRKKQ LDYYQNFLFN WIQKQKELVE QASTN
<b>Source</b>	E.coli
<b>Target Names</b>	MPN_290
<b>Protein Names</b>	Recommended name: Putative type-1 restriction enzyme specificity protein MPN_290 Alternative name(s): S.MpnORFEAP Type I restriction enzyme specificity protein MPN_290 Short name= S protein
<b>Expression Region</b>	1-145
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.