



Recombinant Escherichia coli Probable fimbrial chaperone protein ElfD (elfD)

Product Code	CSB-MP301981ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75856
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	AGKG GIGLAATRLV YSEGEEQISL GVRNTSPDVP YLIQSWVMTP DNKKSADFII TPPLFVLNPA NENLLRIMYI GAPLAKDRET LFFTSVRAVP STTKRKEGNT LKIATQSVIK LFWRPKGLAY PLGEAPAKLR CTSSADMVTV SNPTPYFITL TDLKIGGKVV KNQMISPFDK YQFSLPKGAK NSSVTYRTIN DYGAETPQLN CKS
Source	Mammalian cell
Target Names	elfD
Protein Names	Recommended name: Probable fimbrial chaperone protein ElfD
Expression Region	27-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.