



Recombinant Escherichia coli 1,2-phenylacetyl-CoA epoxidase, subunit A (paaA)

Product Code	CSB-EP301996ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76077
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MTQEERFEQR IAQETAIEPQ DWMPDAYRKT LIRQIGQHAH SEIVGMLPEG NWITRAPTLR RKAILLAKVQ DEAGHGLYLY SAAETLGCAR EDIYQKMLDG RMKYSSIFNY PTLSWADIGV IGWLVDGAAI VNQVALCRTS YGPYARAMVK ICKEESFHQR QGFEACMALA QGSEAQKQML QDAINRFWWP ALMMFGPNDD NSPNSARSLT WKIKRFTNDE LRQRFVDNTV PQVEMLGMTV PDPDLHFDTE SGHYRFGEID WQEFNEVING RGICNQERLD AKRKAWEEGT WVREAALAHA QKQHARKVA
Source	E.coli
Target Names	paaA
Protein Names	Recommended name: 1,2-phenylacetyl-CoA epoxidase, subunit A EC= 1.14.13.n9 Alternative name(s): 1,2-phenylacetyl-CoA epoxidase, catalytic subunit alpha 1,2-phenylacetyl-CoA monooxygenase, subunit A
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.