



# Recombinant Danio rerio T-box transcription factor TBX6L (tbx6l)

<b>Product Code</b>	CSB-BP302106DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P79742
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNHLANNYGY YPQDCRTQYS RMNSAEAEALT SLPVHVSLQD RELWDFSSI GTEMLITKSG RRMFPSCKVT VTGLNPKVKY VVIMDMVPFD NHKYKWNKDC WEVNGSSDPH LPNRFFIHPD SPAPGQKWMQ YPISFHKLKL TNNTLNSNGL VVLHSMHKYQ PRLHIVQSPD PCTPHNPGAY LRFTFPEAAF IAVTAYQNQE ITKLKIDNNP FAKGFRDNGL NRKRFRDKGT QEMQDTRQV KDLDTANECA AGMSQMVEDV DVSVSSSVDC RDTQNSSSVS LNPFISAFTN PSSAGGAAAH QTHTLLSLSN RHFSSPRESN LNSVCAALPV SQLSTGHTSF SRLNPQETHH NSRPKIQLQP PPHSLQCHDL DVRLPLPPKL SRVQLSESAL RNLEMSPLSD CANPRPLTNI LNRSCFRASR PSGKLLPNPP QPEQFLRGSE REIYPAVQEY TDQQFTLNSQ TEHRPHMRPL TEY
<b>Source</b>	Baculovirus
<b>Target Names</b>	tbx6l
<b>Protein Names</b>	Recommended name: T-box transcription factor TBX6L Short name= T-box protein 6L Alternative name(s): T-box transcription factor TBX6 Short name= T-box protein 6
<b>Expression Region</b>	1-473
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.