



Recombinant Danio rerio T-box transcription factor TBX6L (tbx6l)

Product Code	CSB-EP302106DIL-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79742
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MNHLANNYGY YPQDCRTQYS RMNSAEAEALT SLPVHVSLQD RELWDFSSI GTEMLITKSG RRMFPSCKVT VTGLNPKVKY VVIMDMVPFD NHKYKWNKDC WEVNGSSDPH LPNRFFIHPD SPAPGQKWMQ YPISFHKLKL TNNTLNSNGL VVLHSMHKYQ PRLHIVQSPD PCTPHNPGAY LRFTFPEAAF IAVTAYQNQE ITKLKIDNNP FAKGFRDNGL NRKRFRDKGT QEMQDTRQV KDLDTANECA AGMSQMVEDV DVSVSSSVDC RDTQNSSSVS LNPFISAFTN PSSAGGAAAH QTHTLLSLSN RHFSSPRESN LNSVCAALPV SQLSTGHTSF SRLNPQETHH NSRPKIQLQP PPHSLQCHDL DVRLPLPPKL SRVQLSESAL RNLEMSPLSD CANPRPLTNI LNRSCFRASR PSGKLLPNPP QPEQFLRGSE REIYPAVQEY TDQQFTLNSQ TEHRPHMRPL TEY
Source	E.coli
Target Names	tbx6l
Protein Names	Recommended name: T-box transcription factor TBX6L Short name= T-box protein 6L Alternative name(s): T-box transcription factor TBX6 Short name= T-box protein 6
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.